

# KeyTec® TR-FRET STK S2 kinase kit



**CAT.&Size:** A1080004S (1,000 tests)  
A1080004L (20,000 tests)  
A1080004B (100,000tests)

VKEYBIO-01-2025

**Storage at:** 2-8 °C

**For Research Use Only**

**Not For Diagnostic Or Therapeutic Use**

## KeyTec® TR-FRET Serine/Threonine Kinase Kit - S2 Substrate Technical Manual

### 1. Introduction


**KeyTec® TR-FRET STK S2 kinase kit** is designed for measuring the Serine/threonine kinase (STKs) activity. Serine/threonine kinase family comprises multiple members. This assay kit has been validated for activity measurement of over 30 Serine/threonine kinases (For more details, please visit our official website or contact our technical support team). The assay provides a biotinylated substrate, Streptavidin-HX<sup>\*1</sup>, a phospho STK-specific antibody labeled with KeyTec® TR-FRET Solar Eu<sup>\*2</sup> and detection reagent. These components enable the quantification of kinase activity by detecting the phosphorylated substrate.

This assay is based on a sandwich immunoassay method using KeyTec® TR-FRET technology, offering a simple, rapid, highly specific and sensitive, as well as reproducible detection process. The principle, outlined in Figure 1, involves two main phases:

**Kinase Reaction:** The specific substrate is phosphorylated by the target kinase in the optimized reaction buffer with ATP.

**TR-FRET Detection:** Adding an antibody and Streptavidin-HX containing EDTA (to terminate the kinase reaction) that specifically binds to the phosphorylated substrate. This binding brings the TR-FRET donor and acceptor into close proximity, enabling resonance energy transfer (RET) upon excitation and generating a TR-FRET signal. The intensity of TR-FRET signal is directly proportional to the level of phosphorylated substrate.

- \*1 KeyTec® TR-FRET HX: TR-FRET Acceptor
- \*2 KeyTec® TR-FRET Solar Eu: TR-FRET Donor



**Figure1.** The principle of KeyTec® TR-FRET kinase activity detection

## 2. Components

| Components                                     | Storage | A1080004S<br>(1,000 tests <sup>*3</sup> ) | A1080004L<br>(20,000 tests <sup>*3</sup> ) | A1080004B<br>(100,000 tests <sup>*3</sup> ) |
|--|---------|---|--|---|
| <b>p-STK mAb - Solar Eu</b><br>Lyophilized     | 2-8 °C  | 1 vial<br>1,000 tests/vial                | 1 vial<br>20,000 tests/vial                | 5 vials<br>20,000 tests/vial                |
| <b>Streptavidin - HX</b><br>Lyophilized        | 2-8 °C  | 2 vials<br>840 pmoles/vial<br>A1020024S   | 3 vials<br>8.4 nmoles/vial<br>A1020024L    | 4 vials<br>40 nmoles/vial<br>A1020024B      |
| <b>STK Substrate 2 - biotin</b><br>Lyophilized | 2-8 °C  | 1 vial<br>25 nmoles/vial<br>A1080009S     | 1 vial<br>250 nmoles/vial<br>A1080009L     | 3 vials<br>250 nmoles/vial<br>A1080009L     |
| <b>DTT</b><br>Lyophilized                      | 2-8 °C  | 1 vial<br>120 µmoles/vial<br>A1010022S    | 1 vial<br>600 µmoles/vial<br>A1010022L     | 1 vial<br>3.0 mmoles/vial<br>A1010022B      |
| <b>MgCl<sub>2</sub> (1 M)</b>                  | 2-8 °C  | 1 vial<br>1.5 mL/vial<br>A1010024S        | 1 bottle<br>6.0 mL/bottle<br>A1010024L     | 1 bottle<br>30 mL/bottle<br>A1010024B       |
| <b>CaCl<sub>2</sub> (1 M)</b>                  | 2-8 °C  | 1 vial<br>1.5 mL/vial<br>A1010043S        | 1 bottle<br>6.0 mL/bottle<br>A1010043L     | 1 bottle<br>30 mL/bottle<br>A1010043B       |
| <b>Kinase Enzymatic<br/>buffer (5X)</b>        | 2-8 °C  | 1 bottle<br>10 mL/bottle<br>A1010018S     | 1 bottle<br>50 mL/bottle<br>A1010018L      | 1 bottle<br>250 mL/bottle<br>A1010018B      |
| <b>Kinase Detection<br/>buffer (with EDTA)</b> | 2-8 °C  | 1 bottle<br>20 mL/bottle<br>A1010019S     | 1 bottle<br>250 mL/bottle<br>A1010019L     | 2 bottles<br>550 mL/bottle<br>A1010019B     |

<sup>\*3</sup> Tests refer to the number of assay wells that can be performed in 96-well or 384-well plates with 20 µL in total reaction volume. The reagents of the kit are suggested to use as recommended.

### 3. Storage

- ◆ Store all reagents according to the recommended conditions. The products are stable for one year from the date of receipt.
- ◆ After thawing, aliquot the stock into single-use volumes (recommended minimum: 10  $\mu$ L) to avoid repeated freeze-thaw cycles. Store these aliquots at -60 °C and below.

### 4. Required Components (Not Supplied)

| Material                              | Brand          | Catalog           |
|---------------------------------------|----------------|-------------------|
| ATP                                   | Bidepharmatech | BD112724          |
| AMP (5 mM)                            | Sigma-Aldrich  | A1752             |
| cGMP (1 mM)                           | Sigma-Aldrich  | G6129             |
| Lipid activator                       | Sigma-Aldrich  | 20-133            |
| Microplate                            | VKEY-BIO       | M2000102N         |
| Top sealing film                      | VKEY-BIO       | M1000102N         |
| Microplate Reader with TR-FRET module | TECAN          | Infinite® 200 PRO |

## 5. Procedure

### 5.1 Reaction system

| Phase                              | Components                                    | Volume <sup>*4</sup> ( 20 $\mu$ L) |
|------------------------------------|---|------------------------------------|
| Kinase reaction<br>( 10 $\mu$ L)   | Test samples or positive control              | 2 $\mu$ L <sup>*5</sup>            |
|                                    | Kinase  | 4 $\mu$ L <sup>*5</sup>            |
|                                    | Serine/threonine kinase<br>substrate - biotin | 4 $\mu$ L (Pre-mix) <sup>*5</sup>  |
|                                    | ATP   |                                    |
| TR-FRET detection<br>( 10 $\mu$ L) | p-STK mAb - Solar Eu <sup>*6</sup>            | 5 $\mu$ L                          |
|                                    | Streptavidin - HX <sup>*6</sup>               | 5 $\mu$ L                          |

<sup>\*4</sup> Recommended Format: Shallow-well 384-well microplate; For 96-well or 1536-well microplates, proportionally scale the reaction system.

<sup>\*5</sup> The proportions of test sample, kinase, substrate and ATP may be adjusted, as long as the concentrations stay the same.

<sup>\*6</sup> It is recommended to premix p-STK mAb - Solar Eu and Streptavidin - HX , add 10  $\mu$ L pre-mix solution in TR-FRET detection.

### 5.2 Reagent preparation

- ◆ Before use, equilibrate all reagents to RT.
- ◆ **For lyophilized powder, centrifuge to collect the powder to the bottom (850 $\times$ g, 1-2minutes).** It is normal for the fragmentation of lyophilized powder.
- ◆ Use the provided/recommended buffers to prepare sample and detection reagents.
- ◆ Prepare sample and detection reagent according to the kit technical manual.
- ◆ Prepare all reagents immediately before use, unless otherwise specified in the "Working Solution Preparation" section.
- ◆ To avoid the degradation of active components, keep 1X Kinase Enzymatic buffer (containing supplemental components) on ice.
- ◆ To minimize the aggregation or non-specific binding, It is recommended to add surfactants in Kinase Enzymatic buffer (5X). Notes: Kinase Enzymatic buffer (5X)

contained 0.05% BSA.

- ◆ Avoid vigorous mixing of all reagents.

## 6. Working Solution Preparation

### 6.1 Enzyme reaction buffer

- ◆ Prepare 1x enzyme reaction buffer: dilute 1 volume of 5X Kinase Enzymatic buffer and kinase supplements required (e.g.,  $MgCl_2$ ,  $MnCl_2$ , DTT, etc.) with 4 volume of ultrapure water.
- ◆ **DTT stock solution:** Before reconstitution, centrifuge DTT to pellet the powder to the bottom ( $850 \times g$ , 1-2 minutes). Add the recommended volume of ultrapure water to the powder and mix gently as below. Aliquot DTT stock solution into single-use tubes to avoid repeated freeze-thaw cycles. These aliquots can be stored at  $-60^\circ C$ .

| DTT Powder          | Buffer           | Volume           | Stock Concentration |
|---------------------|------------------|------------------|---------------------|
| 120 $\mu$ mole/vial | Ultra-pure water | 120 $\mu$ L/vial | 1 M                 |
| 600 $\mu$ mole/vial | Ultra-pure water | 600 $\mu$ L/vial | 1 M                 |

### 6.2 Test sample or positive control preparation

- ◆ Dilute the test samples to working concentration (5X). For example, if the final test sample concentration is 10 nM in 10  $\mu$ L of kinase reaction system, dilute test samples to working concentration (50 nM) using 1X enzyme reaction buffer and add 2  $\mu$ L/well of test samples. The percentage of DMSO must not exceed 2% for serial dilution.

### 6.3 Substrate and ATP pre-mix solution

- ◆ **ATP stock solution:** Prepare 1 mM ATP stock solution with ultrapure water, aliquot it into single-use tubes to avoid repeated freeze-thaw cycles. These aliquots can be stored at  $-80^\circ C$  for up to 3 months.
- ◆ **Substrate stock solution:** Before resuspension, centrifuge STK substrate 2 - biotin

to pellet the powder to the bottom (850 × g, 1-2 minutes). Add the recommended volume of ultrapure water to the powder and mix gently as below.

| Substrate Powder | Buffer           | Volume      | Stock Concentration |
|------------------|------------------|-------------|---------------------|
| 25 nmoles/vial   | Ultra-pure water | 500 µL/vial | 50 µM               |
| 250 nmoles/vial  | Ultra-pure water | 500 µL/vial | 500 µM              |

- ◆ **Substrate and ATP pre-mix solution(2.5X)** : Prepare 1 mL substrate and ATP pre-mix Solution(2.5X), dilute 50 µL of STK substrate 2 - biotin stock solution and 12.5 µL of ATP stock solution with 937.5 µL of 1x enzyme reaction buffer.

| Component          | Final Conc. µM | Working Conc. µM | Stock Conc. µM | Dilution | Stock Vol. µL | 1X Enzyme Reaction Buffer µL | Total Vol. µL |
|--------------------|----------------|------------------|----------------|----------|---------------|------------------------------|---------------|
| Substrate - biotin | 1.0            | 2.5              | 50             | 20       | 50            | 937.5                        | 1,000         |
| ATP                | 5              | 12.5             | 1,000          | 80       | 12.5          |                              |               |

## 6.4 Kinase solution

- ◆ Dilute the kinase to working concentration (2.5X). For example, if the final kinase concentration is 1 ng/µL in 10 µL of kinase reaction system, dilute 100 ng/µL to working concentration (2.5 ng/µL) using 1X enzyme reaction buffer and add 4 µL/well of kinase solution. **Prepare it immediately prior to addition.**

## 6.5 Detection solution

- ◆ **p-STK mAb - Solar Eu stock solution**: Before reconstitution, centrifuge p-STK mAb - Solar Eu to pellet the powder to the bottom (850 × g, 1-2 minutes). Add the recommended volume of Kinase detection buffer to the powder and mix gently as below.

| Kit Cat.    | Size              | Buffer                  | Volume      | Stock Concentration |
|-------------|-------------------|-------------------------|-------------|---------------------|
| A1080004S   | 1,000 tests/vial  | Kinase Detection buffer | 0.5 mL/vial | 10X                 |
| A1080004L/B | 20,000 tests/vial | Kinase Detection buffer | 1.0 mL/vial | 100X                |

- ♦ **p-STK mAb - Solar Eu working solution(1X)**: Dilute 10X p-STK mAb - Solar Eu stock solution to working concentration(1X). For example, prepare 1 mL working solution(1X), dilute 100  $\mu$ L of p-Tyr mAb - Solar Eu stock solution with 900  $\mu$ L of kinase detection buffer. Working solution(1X) can be stored at 2-8°C for 24 hours or -60°C for 3 months .
- ♦ **Streptavidin - HX stock solution**: Before reconstitution, centrifuge Streptavidin - HX to pellet the powder to the bottom (850  $\times$  g, 1-2 minutes). Add the recommended volume of Ultra-pure water to the powder and mix gently as below.

| Component Cat. | Size            | Buffer           | Volume           | Stock Concentration |
|----------------|-----------------|------------------|------------------|---------------------|
| A1020024S      | 840 pmoles/vial | Ultra-pure water | 42 $\mu$ L/vial  | 20 $\mu$ M          |
| A1020024L      | 8.4 nmoles/vial | Ultra-pure water | 420 $\mu$ L/vial | 20 $\mu$ M          |
| A1020024B      | 40 nmoles/vial  | Ultra-pure water | 2.0 mL/vial      | 20 $\mu$ M          |

- ♦ **Streptavidin - HX working solution(1X, 250nM)**: The molar ratio of substrate-biotin: Streptavidin-HX has been optimized at 8:1. Dilute Streptavidin-HX stock solution to working concentration (1X). For example, prepare 0.8 mL working solution(1X), dilute 10  $\mu$ L of Streptavidin - HX stock solution with 790  $\mu$ L of kinase detection buffer. Working solution(1X) can be stored at 2-8°C for 24 hours or -60°C for 3 months .

- ♦ **(Optional) p-STK mAb - solar Eu and streptavidin - HX pre-mix solution(0.5X):**  
The pre-mix solution should be combined in a 1:1 volume ratio of p-STK mAb-Solar Eu working solution and streptavidin-HX working solution.

| Component                              | Volume                           | Kinase reaction Conc. (10 $\mu$ L in total) | Final Conc. (20 $\mu$ L in total) | Working Conc.        |
|--|----------------------------------|---|-----------------------------------|----------------------|
| Test samples or positive control       | 2 $\mu$ L/well                   | \   | \                                 | \                    |
| Kinase                                 | 4 $\mu$ L/well                   | 0.1 ng/ $\mu$ L                             | \                                 | 2.5 ng/ $\mu$ L      |
| Substrate - biotin                     | 4 $\mu$ L/well (Premix with ATP) | 1000 nM                                     | 500 nM                            | 2500 nM              |
| Streptavidin - HX (1/8 <sup>*7</sup> ) | 5 $\mu$ L/well                   | \   | 62.5 nM <sup>*7</sup>             | 250 nM <sup>*7</sup> |
| p-Tyr mAb - Solar Eu                   | 5 $\mu$ L/well                   | \   | \                                 | 1X                   |

<sup>\*7</sup> Streptavidin - HX / substrate - biotin = 1/8 (Final Conc.) .

## 7. Procedure

- ♦ Follow the steps in the table below.

|                   |        | Kinase Group  | Negative Control   |
|-------------------|--------|---|--|
| Kinase reaction   | Step 1 | 2 $\mu$ L Test sample (5X)<br>4 $\mu$ L Kinase (2.5X)<br>4 $\mu$ L STK substrate 2 - biotin and ATP pre-mix solution (2.5X)   | 2 $\mu$ L Test sample (5X)<br>4 $\mu$ L Enzymatic buffer<br>4 $\mu$ L STK substrate 2 - biotin and ATP pre-mix solution (2.5X) |
|                   | Step 2 | Seal the microplate to prevent liquid evaporation, Incubation (RT, 25 $^{\circ}$ C) for 60 minutes  |  |
| TR-FRET detection | Step 3 | Remove the top sealing film, add detection solution:<br>5 $\mu$ L p- STK mAb - Solar Eu (1X)<br>5 $\mu$ L Streptavidin - HX (1X)<br>or 10 $\mu$ L p-STK mAb - solar Eu and streptavidin - HX pre-mix solution(0.5X) |  |
|                   | Step 4 | Seal the microplate to prevent liquid evaporation, Incubation (RT, 25 $^{\circ}$ C) for 60 minutes or over-night  |  |
|                   | Step 5 | Record the data in microplate reader with top sealing film  |  |



## 8. Data Analysis


- ♦ Calculate the 665 nm/620 nm Ratio (TR-FRET Ratio) and the percentage coefficient of variation (CV %) for each well.

$$\text{TR-FRET Ratio} = \frac{\text{Signal 665 nm}}{\text{Signal 620 nm}} \times 10,000$$

## 9. Summary


### 9.1 Results

|                              |  |
|------------------------------|--|
| <b>PAK4</b>                  | UniProt ID: O96013   |
| <b>Enzyme concentration</b>  | 0.1 ng/μL [ Final concentration in 10 μL buffer ; 1.0 ng/well ]        |
| <b>Incubation time at RT</b> | 60 min   |
| <b>Substrate</b>             | STK Substrate 2 - biotin, 1 μM [ Final concentration in 10 μL buffer ] |
| <b>Streptavidin - HX</b>     | 62.5 nM [ SA/Substrate = 1/8 for final concentration in 20 μL buffer ] |
| <b>Detection antibody</b>    | p-STK mAb - Solar Eu, 1X [ Working concentration ]                     |
| <b>ATP concentration</b>     | 100 μM   |
| <b>Buffer additives</b>      | 5 mM MgCl <sub>2</sub> , 1 mM DTT                                      |



## 9.2 Performance

- The comparative data presented below show the performance of the VKEY-BIO KIT and Company C KIT under the same experimental conditions.



**Note:** Exemplary data shown. Results are instrument-dependent.

## 10. Instrument Model and Setting

|                   |   |
|-------------------|---|
| Vendor            | TECAN   |
| Instrument model  | Infinite® 200 PRO [ Ref. 30050303 ]                     |
| Mode              | Fluorescence Top Reading                                |
| Excitation filter | 320 (25) nm [ Ref. 30094454 ]                           |
| Emission filter 1 | 665 (8.5) nm [ Ref. 30094518 ]                          |
| Emission filter 2 | 620 (10) nm [ Ref. 30094505 ]                           |
| Mirror            | Dichroic 510  |
| Lag time          | 150 µs  |
| Integration Time  | 500 µs  |
| Number of reads   | 5 or user-defined                                       |
| Gain              | 150 or optimal  |
| Z -focus (mm)     | Can be calculated on the well giving the highest signal |