

KeyTec® 384-Well Black Flat Microplates

PS, Solid, Non-treated, No lid



CAT. & Size M2000402N (40 Pcs/Box)

VKEYBIO-01-2024

For Research Use Only

Storage at RT

Not For Diagnostic Or Therapeutic Use

KeyTec® 384-Well Black Flat Microplates


PS, Solid, Non-treated, No lid

Instruction Manual

1. Introduction

KeyTec® 384-Well Black Flat Microplate is designed for Fluorescence, Fluorescence Polarization, Luminescent assays. With square wells designed for standard assay volumes (suitable for 20 - 100 μ L), it ensures easier operation and reduces errors.

2. Dimension



3. Specifications

CAT.	Description
Automation Compatible	Yes
Coating Treatment	Untreated
Color	Black
Detection Method	Fluorescence, Fluorescence Polarization, Luminescent
Material	Polystyrene
Size	40 Pcs/Box
Wells Number	384 well plate
Well Shape	Square well
Working Volume	20 μ L - 100 μ L
Total Volume	120 μ L